

ElectroSim-LSV-RRDE-Analyzer

■ **Live app:** <https://lsv-analysis-ir-compensation-tafel-slope.streamlit.app/>

A Streamlit app for electrochemistry workflows, run from uploaded Excel/CSV files — no coding required:

1. **iR compensation** — correct linear-sweep voltammetry (LSV) data for ohmic drop, using the uncompensated resistance **R_u** fitted from an EIS Nyquist arc.
2. **LSV analysis** — onset potential, overpotential at benchmark current densities (e.g. $j = 10 \text{ mA/cm}^2$), and the Tafel slope of a polarization curve, with RHE reference conversion and mechanistic interpretation (HER, HOR, OER, ORR, CO₂RR, N₂RR, NO₂RR, and more).
3. **Koutecky-Levich (K-L) analysis** — the classic multi-rotation-rate RDE fit for the kinetic current density and electron-transfer number.
4. **ORR / RRDE analysis** — onset/ $E^{1/2}$ /Tafel/electron-number/peroxide-yield from ring-disk data at one rotation rate, plus a multi-rotation-rate ring/disk comparison.

Quick start

```
pip install -r requirements.txt
streamlit run app.py
```

Open the URL shown in the terminal (default `http://localhost:8501`). In the sidebar, pick **Use bundled sample** to try the EIS/LSV tabs immediately with the example workbook in `sample-data/`.

TIFF downloads need a headless Chrome (via `kaleido`); if it's missing, run `plotly_get_chrome`, or use the HTML/CSV download options instead.

How to use it

1 · EIS / Ru Analysis and 2 · LSV iR Correction

1. **Upload** an Excel workbook (Sheet 1 = EIS: Z' , Z'' ; Sheet 2 = LSV: Potential, Current) or two CSV files. Several samples can sit side-by-side in one sheet as repeated column pairs; pick the active one from the sidebar.
2. **Confirm units.** Current unit and R_u unit are auto-detected from the column headers. If one is per-area (e.g. mA/cm^2) and the other isn't, enter the **electrode area (cm^2)** to reconcile them.
3. **EIS / Ru Analysis tab** — R_u is read off the high-frequency intercept of a circle fit to the Nyquist semicircle (adjust the fit range, or enter R_u manually).
4. **LSV iR Correction tab** — pick one or more compensation factors (5–100 %; **85 % recommended**):

```
text E_corrected = E_measured - (factor / 100) · I · Ru
```

View the corrected curve against the raw one, then download the result as CSV. An over-compensation ("fold-back") guard flags any factor that's unsafe and suggests the highest safe one.

3 · LSV Analysis

Independent of the tabs above — its own uploader, units, and reference electrode.

1. **Upload** one or more Potential/Current files (Excel or CSV), one per sample; combine them into a single overlay plot.
2. **Convert to RHE** if your data isn't already on that scale: pick the reference electrode used (SCE, Ag/AgCl variants, Hg/HgO, ...) and the electrolyte pH.
3. Read off the **onset potential** (where |current| first departs from the flat pre-onset baseline) and the **overpotential at one or more benchmark current densities** (e.g. η at $j = 10 \text{ mA/cm}^2$, the standard OER/HER activity benchmark; $j = 2 \text{ mA/cm}^2$ is also common) — reported as $\eta = |E - E_{\text{eq}}|$ for reactions with a known equilibrium potential (HER/HOR/OER/ORR/NO $\blacksquare$ RR/N $\blacksquare$ RR/CO $\blacksquare$ RR), or as the raw potential otherwise.
4. The **linear Tafel region is auto-detected** from the reaction onset; fine-tune it per sample by dragging the slider or box-selecting the region directly on the plot.
5. **Tag each sample's reaction** (HER, HOR, OER, ORR, CO $\blacksquare$ RR, N $\blacksquare$ RR, NO $\blacksquare$ RR, and more) — samples of different reactions can share one overlay, split into per-reaction legend groups.
6. **Group repeat scans of one sample** by giving them the same **Replicate group** name (auto-filled when you upload the same file more than once) — every fitted value (Tafel slope, onset, $\eta@j$, ...) is then also reported as a **mean $\pm$ SD across the group** in its own results section.
7. Read off the **Tafel slope, R^2 , and nearest mechanistic benchmark**, then download the plots (TIFF/HTML) and plotted data (CSV) — publication-styled with Arial fonts and a selectable font size.

4 · K-L (Koutecky–Levich) Analysis

Independent of the tabs above. Needs an RDE polarization curve at **three or more rotation rates** for the same sample (working/disk electrode current only — no ring needed).

1. **Upload** one file per rotation rate (Excel or CSV, just Potential, Current), or a ZIP of a whole data folder for batch loading (see [Loading RRDE/RDE data](#) below).
2. **Convert to RHE** and set the **electrode area** the same way as the other tabs.
3. Pick the **electrolyte's O $\blacksquare$ transport parameters** (diffusion coefficient, kinematic viscosity, bulk O $\blacksquare$ concentration) — presets for 0.1 M KOH, 0.1 M HClO $\blacksquare$, and 0.5 M H $\blacksquare$ SO $\blacksquare$ are built in, or enter your own.
4. At each of several **analysis potentials** (evenly spaced across the range every rotation rate shares), the app fits

$$\text{text } 1/j = 1/j_k + 1/(B \cdot \omega^{0.5}), \quad \omega = 2\pi \cdot \text{rpm}/60, \quad B = 0.62 \, n \, F \, D^{2/3} \, \nu^{(-1/6)} \, C$$

(Koutecký & Levich, *Zh. Fiz. Khim.* **1958**, 32, 1565; Bard & Faulkner, *Electrochemical Methods*, 2nd ed., Wiley, 2001, Ch. 9) — the **intercept** gives the kinetic current density j_k , and the **slope** gives the electron-transfer number n via the Levich constant B . 5. Download the RDE-curve and K-L plots (TIFF/HTML) and the results table (TIFF/CSV).

5 · ORR / RRDE Analysis

Independent of the tabs above. Works with disk-only (RDE) or disk+ring (RRDE) data, one or more samples, one or more rotation rates each.

1. **Upload** data for one or more samples (see [Loading RRDE/RDE data](#) below) and pick which samples to compare.
2. **Convert to RHE**, set the **electrode area**, and — if ring current is available — the **ring collection efficiency N**.
3. Pick a **primary rotation rate** (conventionally 1600 rpm); each sample uses its own closest available rate if it doesn't have that exact one.
4. At that rotation rate, read off each sample's **onset potential**, **half-wave potential** $E_{1/2}$ (the steepest point of the disk curve), and — after removing the mass-transport contribution — the **Tafel slope** (each sample gets its own adjustable fit-range slider).
5. If ring current is available: **electron-transfer number n** and **peroxide yield %H₂O₂**, computed directly from disk + ring current (no multi-rotation-rate fit needed):

```
text n = 4|I_d| / (|I_d| + |I_r|/N) %H2O2 = 200(|I_r|/N) / (|I_d| + |I_r|/N)
```

(Bard & Faulkner, *Electrochemical Methods*, 2nd ed., Wiley, 2001.) 6. If a sample has **several rotation rates**, see them compared in one merged plot — ring current above zero, disk current below, sharing one potential axis (the usual published-RRDE-figure style). 7. Download every plot (TIFF/HTML) and the results/plotted data (TIFF/CSV).

Loading RRDE/RDE data

The K-L and ORR/RRDE tabs share the same loader, in two styles (usable together):

- **Per-sample uploaders** — for each sample, upload either raw per-electrode files (one file per rotation rate *and* electrode, just Potential, Current — the layout most RDE/RRDE instrument software exports, e.g. Disk Current vs Disk Potential (1600 RPM).csv), or a compiled workbook with Potential, Disk current, Potential, Ring current already paired for one rotation rate. Rotation rate and disk/ring role are guessed from the filename and can be corrected by hand.
- **Batch ZIP upload** — zip a whole data folder (one subfolder per sample) and upload it in one go; role and rotation rate are read from filenames automatically, with no per-file correction UI. Use the per-sample uploaders instead for anything that needs fixing by hand.

Good to know

- **LSV data.** Linear-sweep voltammetry records current as the electrode potential is swept linearly — every tab starts from a Potential, Current (or current-density) pair read from your file.
- **Potential → RHE conversion.** If your data is referenced to another electrode (SCE, Ag/AgCl, ...), it's converted with $E(\text{RHE}) = E(\text{measured}) + E^\circ(\text{ref vs NHE}) + 0.0592 \text{ V} \cdot \text{pH}^{-1} \times \text{pH}$. Skip this if your data is already vs RHE.
- **Current & current-density conversion.** Switching the current-unit dropdown (A ↔ mA ↔ μA ↔ nA, or their /cm² counterparts) instantly rescales the values within that family. Converting *between* an absolute current and a density needs the **electrode area**: $\text{current density [A/cm}^2\text{]} = \text{current [A]} / \text{area [cm}^2\text{]}$. The same area also reconciles Ru between Ω and Ω·cm² so $I \cdot R_u$ always resolves to volts.

- **Privacy.** Uploaded files are processed in memory only — nothing is written to disk — and a **Clear / reset** button wipes the session. An optional password gate can be set via a `app_password` secret.
- Full technical notes (data format, security hardening, using the Python API without the GUI) are in the module docstrings under [scripts/modules/](#).

Citation

If you use this app (including the public Streamlit deployment) in your work, please cite it. Machine-readable metadata is in [CITATION.cff](#); a plain-text form:

Kumar, R. (2026). *ElectroSim-LSV-RRDE-Analyzer* (v1.1.0) [Computer software]. North Carolina Central University. <https://github.com/rajeev4187/LSV-Analysis-iR-compensation-Tafel-slope>

```
@software{kumar_electrosim_lsv_rrde_2026,  
  author = {Kumar, Rajeev},  
  title  = {ElectroSim-LSV-RRDE-Analyzer},  
  version = {1.1.0},  
  year   = {2026},  
  url    = {https://github.com/rajeev4187/LSV-Analysis-iR-compensation-Tafel-slope}  
}
```

License

Released under the MIT License — see [LICENSE](#).